

ELECTRIC two Section Massage Table



massage-table-electric.jpg



ELECTRIC two Section Massage Table



ELECTRIC two Section Massage Table

?????????: ??? ?????????? ????
????
1225,12 €

1225,12 €

✓ Preorder - προπαραγγελία

[???? ??? ??????? ??? ?? ?????](#)
?????????: [Quirumed](#)

????????

New, modern and versatile two Section Massage Table with electric lifting system, moveable backrest, face cap and thick padding. Professionals, due to its easy functioning and features, will feel in total control and will work in an agile and comfortable way. At the same time, customers will feel totally safe and comfortable.

CHARACTERISTICS:

- Structure reinforced and made in painted steel tube.
- Height adjustment by electric motor, manual remote included.
- Backrest, with padded hole plug, adjustable in inclination from 0 to 50° by hydraulic system.
- Upholstered in fireproof PU leather.
- Maximum weight supported in work up to 200 Kg approx.
- High density padding.
- Roll holder.
- Available Colours :White, Blue , Silver , Black , Green , Came , Red

TECHNICAL DATA:

- **Dimensions:**
 - Width: 70 cm
 - Length: 188 cm
 - Height adjustable: 52-86 cm
- **Padding thickness:** 8 cm
- **Dimensions upper body:**
 - Width: 70 cm
 - Length: 73 cm
- **Dimensions lower body:**
 - Width: 70 cm
 - Length: 115 cm

ELECTRIC two Section Massage Table

- **Dimensions facial hole:**
 - Width: 14 cm
 - Length: 16 cm
 - **Weight:** 55 kg
 - **Packing:** 187 × 74 × 52 cm, weight 75 kg.
- A must-have physiotherapy table for physiotherapists and sports teams

?????? ??? ??????????:1