


## NubisPro – Portable physiotherapy table





NubisPro – Portable physiotherapy table

?????????: ??? ????????????? ?????

????

1395,00 €

1395,00 €

✓ Preorder - προπαραγγελία

[???? ??? ??????? ???? ?? ??????](#)

?????????????: [Nubis](#)

?????????

## Most portable physiotherapy table

NubisPro is a professional therapeutic table that can be folded into a backpack. This is equipped with lined straps that evenly distribute the weight to the body. Because of that NubisPro could be carried hands-free and it is much easier to carry it up and down the stairs. The backpack also features rolling wheels to move faster on the city streets. Due to the compact size, it can fit into the smallest car trunk.

## Quickly prepared

It only takes 2-3 minutes to be completely prepared to help an athlete.

## Extreme weight capacity

Base of the table is made from superior and extremely solid aluminium construction, which gives to NubisPro an amazing weight capacity and can easily hold both, athlete and therapist.

### **Therapist friendly**

Transporting NubisPro is comfortable and simple, even by foot or by bus.

### **Superior material**

High-quality PU surface is oil and waterproof, flame retardant and steady against blood, urine and disinfectants.

### **Easy cleaning**

A NubisPro changeable cover is made of a very durable material, which can be easily cleaned and disinfected.

### **Patient comfort**

Pressure based system, shifting headrest and very flexible and elastic top surface material gives patient ultimate comfort during treatment.

### **Pleasure to use**

NubisPro is a pleasure to work with, thanks to easy system for height adjustment.

### **Attractive**

NubisPro is pleasantly rounded from all sides. There are covers available in several attractive colors for every taste.

Includes:

- Mattress
- Headrest
- Frame and legs

- electric pump
- backpack

**Specs**

Mattress size: 190 x 75 x 20 cm

Made : double wall fabric / PU

Inflation time: 35 s

Legs : Aluminium

Height adjustment : 60-85 cm

Weight capacity: 250 kg

Table weight : : 10,2 kg

Backpack measures : 70 x 45 x 20 cm

Transport weight: 14,7 kg

**Available Colors :**

???, ????????, ?????, ?????? ???, ?????? ???



?????? ? ??:1