



Pagani Cryo-T Shock ,desktop - cryo hot therapy





Pagani Cryo-T Shock ,desktop - cryo hot therapy

?????????: ??? ????????????? ????
???

7999,24 €
7499,52 €
-499,72 €

✓ Preorder - προπαραγγελία

[???? ???? ?????? ??? ?? ?????](#)

?????????: [Elettronica Pagani](#)

????????

Pagani Cryo-T Shock ,desktop - cryo hot therapy

Naturally innovative, ROLAND TOWER CRYO T-SHOCK is a very compact device offering pioneering solutions allowing you to get the best from each treatment. Temperature could be widely adjustable from -5°C to 41°C. TOWER CTYO T-SHOCK is mainly used into the Rehabilitation and Sports Medicine or into the Aesthetics for performing face-lifts.

TOWER CRYO SHOCK allows you to "attack" and treat immediately also recent traumas, on which in most cases not you could intervene and you should ask the patient to be reviewed only after 3-4 days, forcing him to actually endure the pain and a slow regression of the consequences of the trauma. CRYO SHOCK allows you to intervene promptly, to treat pain and reduce blood supply, edema, inflammation.

The patient gets a prompt regression of pain and recovery times are drastically reduced.

Further to various applications into the rehabilitation field, TOWER could be largely suitable for performing all treatments belonging to face-lifts purposes. Applicator size, temperature and patient comfort are absolutely suitable for getting very good performances with care and safety.

For treating patients at the Clinic or at home, TOWER CRYO T-SHOCK is the perfect device for all operators requiring for a compact device, easy to use, convenient in price.

Main features:

- Universal power supply for 110 and 230 V, 50/60 Hz
- Electric Class 1, Type B
- Power consumption: max. 120 VA
- Cold temperature: max – 5°C
- Hot temperature: max. 41°C
- 4" resistive coloured display
- Manual hand-piece with radiating area of 45 mm large
- Highest performances thanks to real-time temperature control with sensor attesting skin's temperature
- Standard accessories: main unit, hand-piece, power cord with "shucko" plug, cryo-gel 1 Kg, operating manual
- On inquiry: adaptor 30 mm large, carrying case-

??????? ??? ??????????:1